



TeamBest™ ASIA

A TEAMBEST GLOBAL COMPANY

www.teambest.com • www.teambest.in

Best™ Integrated Brachytherapy Solutions

You don't go to a shoe store and cut your feet to fit the shoe they have in stock ... then why compromise when it comes to improving the clinical outcome for your patients?



Best™
Universal
Chair/Table


Best™ HDR Remote
Afterloader



BestMick Applicator



BestMick 10/15/20/25 Cartridges



Best™ Kobold Tandem
& Ring Applicator


Best™ Kobold Fletcher
Tandem & Ovoid Applicator

Best™ Kobold Henschke
Tandem & Ovoid Applicator

Best™ Compact
Sonalis™ Ultrasound
Imaging System



Best™ X-ray Marker Cables




mobileMOSFET Wireless
Dosimetry System



Portable
Dosimetry
System

Best™ Templates



Best™
Stepper/Stabilizer



Best™
Treatment
Planning
System



Best™ Esophageal
Brachytherapy Applicator




Best™ Breast Double-Balloon
Brachytherapy Applicator

Best™ Localization Needles with I-125 Seeds



Best™ Fiducial Markers



Best™ Iodine-125 Seed




Best™ Palladium-103 Seed



Stranded

Loose

Best™ Flexi Needles



15G Flexi
Needle

16G Flexi
Needle

17G Flexi
Needle

15G Flexi Needle with
SS Sleeve (Click Fit)

Flexi Needle with
Female Hub



Huestis CVO-2000
Warming Oven for
Thermoplastics



We ship your order within 24 hours, 7 days
a week, anywhere in the world!

Standard or customized marker configuration made
to your exact specifications, sterile or non-sterile!




Best[™] medical international

A TEAMBEST GLOBAL COMPANY

www.teambest.com • www.teambest.in



**Best[™] Raycell
Mk1 Blood
& Research
Irradiator**




**Best[™] Raycell
Mk2 Blood
& Research
Irradiator**



GammaBeam[™]
saving lives

**Best[™] GammaBeam[™] 100/300 Equinox
Teletherapy System with IMRT & IGRT**



UPGRADE KIT for all old Theratron units, 80 or 100 cm including IMRT capabilities w/built in or external MLC



UPGRADE includes:

- Removing all old controls, electronics and installing a new control system and covers
- Replacing the old collimator system with the new Equinox collimator
- Replacing the old treatment table with the new Avanza[™] Table
- Retaining the head rotation capability is optional

UPGRADE features:

- Calculated Arc Speed
- Graphical Control System
- Asymmetric Jaws (optional)
- R&V System Ready (optional)
- Service Log Files
- On-Board Verification
- Motorized Wedge (optional)
- Collision Detection (optional)

**Best[™]
Gammacell[®]
1000/3000 Blood
& Research
Irradiator**






TeamBest™ ASIA

A TEAMBEST GLOBAL COMPANY

www.teambest.com • www.teambest.in




NEW! Best™ BG-95 Sub-Compact Self-Shielded Cyclotron for production of: ^{18}F FDG, Na^{18}F , ^{18}F -MISO, ^{18}F FLT, ^{18}F -Choline, ^{18}F -DOPA, ^{18}F -PSMA, ^{13}N and ^{68}Ga



Best Model B25p Cyclotron

Best Model 6-15 MeV Compact High Current/ Variable Energy Proton Cyclotron

Best Model B35ADP Alpha/ Deuteron/Proton Cyclotron





Best™ 70 MeV Cyclotron at INFN, Legnaro, Italy

B100 Cyclotron	7.5 MeV	<ul style="list-style-type: none"> Capable of producing: ^{18}FFDG and Na^{18}F Single or batch dose production Integrated self-shielded cyclotron, chemistry module and FDG QC module Complete production lab in a 5 x 5 meter area
BG-95 Cyclotron	1-9.5 MeV	<ul style="list-style-type: none"> Low energy, self-shielded compact system capable of producing: ^{18}FFDG, Na^{18}F, ^{18}F-MISO, ^{18}FFLT, ^{18}F-Choline, ^{18}F-DOPA, ^{18}F-PSMA, ^{13}N and ^{68}Ga
Best Cyclotrons	1-3 MeV	Deuterons for materials analysis*
	70-200 MeV	For Proton Therapy*
	3-90 MeV	High current proton beams for neutron production and delivery*
B6-15 Cyclotron	1-15 MeV	Proton only, capable of high current up to 1000 Micro Amps, for medical radioisotopes
B25 Cyclotron	20, 15-25 MeV	Proton only, capable of high current up to 1000 Micro Amps, for medical radioisotopes
B25u-35adp Cyclotron	25-35 MeV	Proton or alpha/deuteron/proton, capable of high current up to 1000 Micro Amps, for medical radioisotopes
B35 Cyclotron	15-35 MeV	Proton only system for medical radioisotopes production
B70/70adp Cyclotron	35-70 MeV	Proton only or alpha/deuteron/proton systems, capable of high current up to 1000 Micro Amps, for medical radioisotopes

Ion Rapid Cycling Medical Synchrotron (iRCMS) 200-400 MeV Variable Energy for Proton to Carbon Heavy Ion Radiation Therapy

Racetrack Synchrotron—Smaller Area Footprint



Best™ Model 200p Cyclotron for Proton Therapy*

Best[™] nomos[®] Compact Sonalis[®] Ultrasound System



The Best[™] Sonalis[®] Ultrasound Imaging System provides superior visualization of HDR, LDR, RF or Cryosurgical procedures. Our patented SimulView[™] Technology provides simultaneous “live” views of the prostate in both planes, thereby increasing treatment accuracy and precision.



Additional Features

- Capable of supporting over 20 different probe configurations, this system can be used in almost every discipline
- 21.5" touch screen technology
- Superior HD image resolution for improved procedure accuracy, speed and physician confidence level
- Patented probe design
- Sagittal array provides for 140 mm length of view encompassing the bladder, prostate and perineal tissue
- Advanced modular software design provides for future upgrade path
- Convertible to desktop or portable without losing any functionality
- Dual battery backup built-in to protect from power loss

*Certain products shown are not available for sale currently.

Sonalis[®] Transducers

We offer probes for the most common types of procedures and imaging, all with the enhanced image quality and high resolution that you would expect from the **BEST** name in ultrasounds!

- 8L2A Linear Array 
- 12L5A Linear Array 
- 14L3 Linear Array 
- 15LW4 Linear Array 
- 15LA Linear Array 
- 15L4A Linear Array 
- 16L5 Linear Array 
- 8V3 Phased Array 
- 4V2A Phased Array 
- 5C2A Curved Array 
- 9MC3 Curved Array 
- 8EC4A Endocavity 
- XY-BI-Plane Phased Array 
- 10EC4A Endocavity 
- 10BP4 Bi-Plane 
- 8BP4 Bi-Plane 
- 8TE3 Trans-esophageal 
- Pedoff 
- 16HL7 High Frequency Linear Array 

TeamBest Global Companies ©2024

Best Nomos • 1 Best Drive • Pittsburgh, PA 15202 USA • tel 412 312 6700 or 800 70 NOMOS • fax 412 312 6701 • www.nomos.com

AccuBoost[®] DIGITAL

PARTIAL BREAST BRACHYTHERAPY SOLUTION

AccuBoost[®] Digital—A digital mammography-guided Non-Invasive Breast Brachytherapy (NIBB) equipment solution. This complete, turn-key solution provides all the components needed for this advanced treatment, bundled together in a streamlined offering.



HUESTIS MEDICAL
A TEAMBEST GLOBAL COMPANY

One Industrial Way | Tyngsboro, MA 01879
1-800-972-9222

www.accuboot.com
www.huestis.com
www.teambest.com