

Best[™] nomos[®] Compact Sonalis[®] Ultrasound System



The Best[™] Sonalis[®] Ultrasound Imaging System provides superior visualization of HDR, LDR, RF or Cryosurgical procedures. Our patented SimuView[™] Technology provides simultaneous “live” views of the prostate in both planes, thereby increasing treatment accuracy and precision.



Additional Features

- Capable of supporting over 20 different probe configurations, this system can be used in almost every discipline
- 21.5" touch screen technology
- Superior HD image resolution for improved procedure accuracy, speed and physician confidence level
- Patented probe design
- Sagittal array provides for 140 mm length of view encompassing the bladder, prostate and perineal tissue
- Advanced modular software design provides for future upgrade path
- Convertible to desktop or portable without losing any functionality
- Dual battery backup built-in to protect from power loss

*Certain products shown are not available for sale currently.

Sonalis[®] Transducers

We offer probes for the most common types of procedures and imaging, all with the enhanced image quality and high resolution that you would expect from the **BEST** name in ultrasounds!

- 8L2A Linear Array 
- 12L5A Linear Array 
- 14L3 Linear Array 
- 15LW4 Linear Array 
- 15LA Linear Array 
- 15L4A Linear Array 
- 16L5 Linear Array 
- 8V3 Phased Array 
- 4V2A Phased Array 
- 5C2A Curved Array 
- 9MC3 Curved Array 
- 8EC4A Endocavity 
- XY-BI-Plane Phased Array 
- 10EC4A Endocavity 
- 10BP4 Bi-Plane 
- 8BP4 Bi-Plane 
- 8TE3 Trans-esophageal 
- Pedoff 
- 16HL7 High Frequency Linear Array 


TeamBest Global Companies ©2024

Best Nomos • 1 Best Drive • Pittsburgh, PA 15202 USA • tel 412 312 6700 or 800 70 NOMOS • fax 412 312 6701 • www.nomos.com

AccuBoost[®] DIGITAL

PARTIAL BREAST BRACHYTHERAPY SOLUTION

AccuBoost[®] Digital—A digital mammography-guided Non-Invasive Breast Brachytherapy (NIBB) equipment solution. This complete, turn-key solution provides all the components needed for this advanced treatment, bundled together in a streamlined offering.



HUESTIS MEDICAL
A TEAMBEST GLOBAL COMPANY

One Industrial Way | Tyngsboro, MA 01879

1-800-972-9222

www.accuboot.com
www.huestis.com
www.teambest.com